



Luer-Lock Connector

11.5⁺³

77±2

1.4305

340±3

1.4301

POM (black)

1.4305

250±3

D-D

B-B

GS1-DataMatrix:
 (01)0426046384xxxx(21)Axxxx
 (01) -> GTIN-Number
 (21) -> Serialnumber
 (xxxx represents actual number)

Sterilisation Tray
 Ref #:172099
 autoclavable, max 134°C
 life time > 500 cycles

for effective reprocessing of:
 Universal Joint, Endoscope Clamp, Tension Sleeves (D4, D5, D10, D12,7),
 Joystick LH, Joystick RH and Mounts (1A, 1B, 2A, 2B)

∇ = ∇	
∇ = ∇	Ra 6.3
∇ = ∇	Ra 3.2
∇ = ∇	Ra 1.6
∇ = ∇	Ra 0.8

Pro/E drawing - do not change		Material:		Part number:		State	
Tolerances DIN ISO 2768-mk		scale:	1:2	-		1308-15000-03	
		weight:	1.343 kg	-			
		Date	Name	Part Name:			
		drawn:	25.03.19	pk			
		chkd:	NYX	NYX			
Surface:				Siebkorb			
Proprietary data, company confidential. All rights reserved.				page:2 of:2			

