


Joystick LH (left hand)
Ref #: 242695

autoclavable, max 134°C
life time max. 150 cycles

connectable with:
Mount 1A 1B
Mount 2A 2B



| | | | | | | | | |
|--|---------|-----------|--|--------------------------------------|--|-------------------------------------|--|--|
| Pro/E drawing - do not change | | Material: | | Part number: | | State | | |
| Tolerances DIN ISO 2768-mk | scale: | 1:1 | | 1828-30000-00 | | | | |
| | weight: | 0.099 kg | | | | | | |
| | Date | Name | | Part Name: Joystick-LH-IIS | | AKTORmed™ ROBOTIC SURGERY | | |
| | drawn: | JL | | | | | | |
| | chkd: | NYX | | | | | | |
| Surface: | | | | | | | | |
| Proprietary data, company confidential. All rights reserved. | | | | | | page:2 of:2 | | |