




GS1-DataMatrix:

(01)0426046384xxxx(21)Axxxx  
 (01) -> GTIN-Number  
 (21) -> Serialnumber  
 (xxxx represents actual number)



Mount 1B  
 Ref #: 172065  
 autoclavable, max 134°C  
 life time > 500 cycles

opening width 2,5 ... 19,5 mm (+/- 1)



✓	=	✓
W	=	Ra 6.3
X	=	Ra 3.2
Y	=	Ra 1.6
Z	=	Ra 0.8

Pro/E drawing - do not change		Material:		Part number:	State	
Tolerances DIN ISO 2768-mk	scale:	2:1		1308-31200-03		
	weight:	0.050 kg				
	Date	Name	Part Name:			
	drawn:	12.12.2023	AZ			Mount 1B
	chkd:	NYX	NYX			
Surface:	Proprietary data, company confidential. All rights reserved.			page:2	of:2	