



GS1-DataMatrix:

(01)0426046384xxxx(21)Axxxx
 (01) -> GTIN-Number
 (21) -> Serialnumber
 (xxxx represents actual number)



Mount 2B
 Ref #: 172068
 autoclavable, max 134°C
 life time > 500 cycles

	=	
W	=	Ra 6.3
X	=	Ra 3.2
Y	=	Ra 1.6
Z	=	Ra 0.8

Pro/E drawing - do not change		Material:		Part number:	State	
Tolerances DIN ISO 2768-mk	scale:	1:1		1308-32100-04		
	weight:	0.028 kg				
	Date	Name	Part Name:			
	drawn:	12.12.2023	AZ			Mount 2B
	chkd:	NYX	NYX			
Surface:	Proprietary data, company confidential. All rights reserved.		page:2 of:2			

Model: 1308-32100_ANBINDUNG-L | Drawing: 1308-32100_ANBINDUNG-L | Size: A4 | 1:1