



Joystick RH (right hand)  
 Ref #:182198  
 autoclavable, max 134°C  
 life time max. 150 cycles

connectable with:  
 Mount 1A 1B  
 Mount 2A 2B

	±0.2	✓ = ✓
	±0.2	✓ = ✓ Ra 6.3
		✓ = ✓ Ra 3.2
		✓ = ✓ Ra 1.6
		✓ = ✓ Ra 0.8

Pro/E drawing - do not change		Material:	Part number:	State	
Tolerances DIN ISO 2768-mk	scale:	1:1	1308-35000-08		
	weight:	0.099 kg			
	Date	Name	Part Name: <b>Joystick-RH</b>		
	drawn:	11.04.24			JL
	chkd:	NYX			NYX
Surface:			<b>AKTORMed™</b> ROBOTIC SURGERY		
Proprietary data, company confidential. All rights reserved.					page:2 of:2