





GS1-DataMatrix:
 (01)0426046384xxxx(21)Axxxx
 (01) -> GTIN-Number
 (21) -> Serialnumber
 (xxxx represents actual number)



✓	=	✓
w	=	✓ Ra 6.3
x	=	✓ Ra 3.2
y	=	✓ Ra 1.6
z	=	✓ Ra 0.8

Tension sleeve D12.7
 Ref #:161773
 autoclavable, max 134°C
 life time > 100 cycles
 suitable for endoscope Ø11.5....12.7
 Color code: green (yellow)

Pro/E drawing - do not change		Material: PPSU MT	Part number: 1206-36107-03	State
Tolerances DIN ISO 2768-mk	scale: 3:1 weight: 0.005 kg			
	Date	Name	Part Name:	
	drawn: 21.03.19	RG	Spannhülse D12,7	
	chkd: NYX	NYX		
Surface:	Proprietary data, company confidential. All rights reserved.		page:2 of:2	