


NYX: managed document! Number: 1515-10002 Index: 01, freigegeben / released: Andreas Lorenz 11.08.2022 11:46:39
 Printed versions of this document are for information only!



Y-Line
 REF#:171895

Cable length 30 +/- 1 cm

Permits simultaneous connection of VOICE CONTROL
 and JOYSTICK to the SOLOASSIST II

	=	
	=	
	=	
	=	

Pro/E drawing - do not change		Material:	Part number:	State	
Tolerances DIN ISO 2768-mk	scale:	2:1	1515-10002-01		
	weight:	0,094 kg			
	Date	Name	Part Name: Y-LINE		
	drawn:	12.12.2023			AZ
	chkd:	NYX			NYX
Surface:					
Proprietary data, company confidential. All rights reserved.					page:2 of:2