



GS1-DataMatrix:

(01)0426046384xxxx(21)Axxxx
 (01) -> GTIN-Number
 (21) -> Serialnumber
 (xxxx represents actual number)



Mount 1A
 Ref #:172066
 autoclavable, max 134°C
 life time > 500 cycles

opening width 2,5 ... 14,5 mm (+/- 1)

	=	
	=	
	=	
	=	
	=	

Pro/E drawing - do not change		Material:		Part number:	State	
Tolerances DIN ISO 2768-mk	scale:	2:1		1308-31000-05		
	weight:	0.043 kg				
	Date	Name	Part Name:			
	drawn:	12.12.2023	AZ			Mount 1A
	chkd:	NYX	NYX			
Surface:	Proprietary data, company confidential. All rights reserved.		page:2 of:2			