



Joystick LH SYK (left hand for Stryker)
Ref #: 212543

autoclavable, max 134°C
life time max. 150 cycles

Power Supply:
24V DC, max. 1,2 W

connectable with:
Mount 1A 1B
Mount 2A 2B

	=	
	=	Ra 6.3
	=	Ra 3.2
	=	Ra 1.6
	=	Ra 0.8

Pro/E drawing - do not change		Material:	Part number:	State
Tolerances DIN ISO 2768-mk	scale:	1:1	1308-33000-00	
	weight:	0.098 kg		
	Date	Name	Part Name: Joystick-LH-SYK	
	drawn:	17.09.2021 AZ		
	chkd:	NYX NYX		
Surface:	Proprietary data, company confidential. All rights reserved.		page:2 of:2	

