





GS1-DataMatrix:

(01)0426046384xxxx(21)Axxxx
 (01) -> GTIN-Number
 (21) -> Serialnumber
 (xxxx represents actual number)



Adapter 2
 for use with Olympus handle
 Ref #: 172088
 autoclavable, max 134°C
 life time min. 100 cycles

connectable with:
 Mount IB

Color: RAL7024 graphit grey

	=	
W	=	Ra 6.3
X	=	Ra 3.2
Y	=	Ra 1.6
Z	=	Ra 0.8

Pro/E drawing - do not change		Material:		Part number:		State	
Tolerances DIN ISO 2768-mk	scale:	2:1		-	1308-20300-04		
	weight:	0.005 kg					
	Date	Name		Part Name:			
	drawn:	AZ					
	chkd:	NYX					
Surface:				Adapter-2			
Proprietary data, company confidential. All rights reserved.				page:2 of:2			

Model: 1308-20300_ADAPTER-2 | Drawing: 1308-20300_ADAPTER-2 | Size: A4 | 2:1